

AFE General Purpose 1 ADC 16bit 3.3V 80-Pin LQFP T/R

Manufacturer:	Texas Instruments, Inc	
Package/Case:	LQFP80	AFE4300PNR Image Images are for reference only
Product Type:	Data Conversion ICs	Inquiry
RoHS:	RoHS Compliant/Lead free W	
Lifecycle:	Active	

General Description

The AFE4300 is a low-cost analog front-end incorporating two separate signal chains: onechain for weight-scale (WS) measurement and the other for body composition measurement (BCM)analysis. A 16-bit, 860-SPS analog-to-digital converter (ADC) is multiplexed between both chains. The weight measurement chain includes an instrumentation amplifier (INA) with the gain set by anexternal resistor, followed by a 6-bit digital-to-analog converter (DAC) for offset correction, and a circuit to drive the external bridge/load cell with a fixed 1.7 V for ratiometric measurements.

The AFE4300 can also measure body composition by applying a sinusoidal current into thebody. The sinusoidal current is generated with an internal pattern generator and a 6-bit, 1-MSPSDAC. A voltage-to-current converter applies this sinusoidal current into the body, between twoterminals. The voltage created across these two terminals as a result of the impedance of the body measured back with a differential amplifier, rectified, and the amplitude is extracted and measured by the 16-bit ADC.

The AFE4300 operates from 2 V to 3.6 V, is specified from 0°C to +70°C, and is available in a LQFP-80 package.

Key Features

Weight-scalefront-end: Supports up to four load cell inputs

On-chip load cell 1.7-Vexcitation voltage for ratiometric measurement

68-nVrms Input-Referred noise(0.1 Hz to 2 Hz)

Best-fit linearity: 0.01% offull-scale

Weight-scale measurement : $540 \mu A$

Body compositionfront-end: Supports up to three tetra-polar complex impedancemeasurements

6-Bit, 1-MSPS sine-wave generation digital-to-analog converter(DAC)

247.5- μ Arms, ±20% Excitation source

0.1-Ω Measurement RMS noise in 2-HzBW

Body composition measurement : 970µA

Analog-to-digital converter(ADC): 16 Bits, 860 SPS

Supply current: 110 µA

All trademarks are the property of their respective owners.

Description

The AFE4300 is a low-cost analog front-end incorporating two separate signal chains: onechain for weight-scale (WS) measurement and the other for body composition measurement (BCM)analysis. A 16-bit, 860-SPS analog-to-digital converter (ADC) is multiplexed between both chains. The weight measurement chain includes an instrumentation amplifier (INA) with the gain set by anexternal resistor, followed by a 6-bit digital-to-analog converter (DAC) for offset correction, and a circuit to drive the external bridge/load cell with a fixed 1.7 V for ratiometric measurements.

The AFE4300 can also measure body composition by applying a sinusoidal current into thebody. The sinusoidal current is generated with an internal pattern generator and a 6-bit, 1-MSPSDAC. A voltage-to-current converter applies this sinusoidal current into the body, between twoterminals. The voltage created across these two terminals as a result of the impedance of the body measured back with a differential amplifier, rectified, and the amplitude is extracted and measured by the 16-bit ADC.

The AFE4300 operates from 2 V to 3.6 V, is specified from 0°C to +70°C, and is availablein a LQFP-80 package.





🟠 AVAQ

🟠 AVAQ

Recommended For You

AFE5807ZCF Texas Instruments, Inc BGA

AFE2124E Texas Instruments, Inc SSOP48

AFE4403YZPT

Texas Instruments, Inc DSBGA36

AFE4400RHAT Texas Instruments, Inc

VQFN40

AFE4404YZPT Texas Instruments, Inc DSBGA15

AFE1205E Texas Instruments, Inc XX

AFE1103E Texas Instruments, Inc SSOP

AFE4403YZPR Texas Instruments, Inc DSBGA36

AFE4490RHAT Texas Instruments, Inc QFN

AFE5808AZCF Texas Instruments, Inc BGA AFE1104E

Texas Instruments, Inc SSOP

AFE5818ZBV

Texas Instruments, Inc BGA

AFE4404YZPR Texas Instruments, Inc DSBGA15

AFE4405YZR

Texas Instruments, Inc DSBGA

AFF5812ZCF Texas Instruments, Inc BGA135